



**MEIJI
TECHNO**

MG900/2X

Equipped with next generation, surface mount LEDs, these innovative designs create the new standard for illuminated magnifiers.

7" x 5.25"
Rectangular
LED Magnilite®



MG900/2XBLK

7" x 5.25"
Rectangular
LED Magnilite®



MG900/2XRED

Lighting Features

- 4x the lumen output of a standard magnifier
- Reduced glare, heat and energy consumption
- Multi-angle LEDs for shadowless lighting at a full range of focal lengths and heights
- User controlled LED groups shine desired level of light exactly where required
- Superb dimming software via convenient, reliable membrane switch
- Built in energy management system & lighting timer

Optical Features

- Precision ground, crown white optical grade glass
- Largest lenses in the industry
- Standard 4 Diopter (2x) lens

Optional Finishes : MG900/2XBLK MG900/2XSIL MG900/2XBLU MG900/2XRED



Made in the U.S.A.

5895 Rue Ferrari San Jose, CA 95138
1 (800) 832-0060 | Fax: (408) 226-0900
www.meijitechno.com