## User-friendly, Affordable

## Binocular Stereo Microscope EM-20 SERIES



MEIJI TECHNO
Robust metal body with smooth rack & pinion

- ·Light weight, compact body, easy to carry with hand grip
- ·Stylish, low-profile light base

Fresh design. Ergonomic focus knobs

- ·LED features a long life with low power consumption
- •Transmitted and incident LED illumination equipped with variable intensity
- Equipped with Rechargeable batteries, cordless
  It enables use of the instrument without an outlet.
- · EM-22 has 1X, 2X, 3X Turret
- · EM-23 has 1X, 2X, 4X Turret
- · Turret Style Objective Changer





MODEL	EM-22	EM-23
Magnification	10X, 20X, 30X	10X, 20X, 40X
Objectives	Turret 1X/2X/3X	Turret 1X/2X4X
Eyepieces	Paired WF10X/20mm	
Viewing head	Fixed Binocular 45 degree inclined	
	Dioptric adjustment on left eyetube, Interpupillary 55-75mm	
Stand	Rack & pinion stand with Hand grip, Charging function with rechargeable batteries	
Focusing block	Viewing head in up-and-down motion	
Working distance	60 mm	
Illumination	LED incident light/transmitted light, intensity control	
Accessories	Paired eyeshields, AC adapter, dust cover	

USA Branch: URL: http://www.meijitechno.com

Meiji Techno America

5895 Rue Ferrari, San Jose, Ca. 95138

Phone: 408-226-3454 Fax: 408-226-0900

Emai: info@meijitechno.com



Manufacturer: URL: http://www.meijitechno.com 322-1, Chikumazawa, Miyoshi-machi, Iruma-gun, Saitama 354-0043, Japan

Phone: +81-49-259-0111 Fax: +81-49-259-0113

Email: meiji@meijitechno.co.jp